

# User manual

# MR902





## ***Table of contents***

Warning and safety instructions:	3
Product description	5
Installation	5
Backup	6
Commissioning	6
Main menu	8
Main menu - Receiver	12
Main menu - Settings	13
Set response time	14
Edit name	14
Size of the system	15
Maintenance	18
Cleaning	18
Spare parts and accessories	18
Technical data	19
Explanation of symbols and approvals	20

## **Warning and safety instructions:**



- The battery in this product complies with the relevant IEC safety standards for batteries.
- Please read the safety instructions on the power supply before use.
- The power supply unit must meet the relevant IEC safety standards.
- The battery must not be replaced, and the product must be sent for repair.
- The product must not be disassembled. Risk of electric shock!
- Read the intended use before putting it into service.
- The product is not water-resistant unless noted in product specification.

## **Contraindications**

In general, the product cannot be used if:

- The disabled client/patient is mentally ill.
- The disabled client/patient is critically ill.
- The client/patient is unable to use the transmitters.

## **Lifetime after mounting of the device**

The battery, if applies, must be replaced when the battery replacement information appears on the display and according to the user manual.

Lifetime is evaluated in relation to the pressure button. This is made according to the supplier's specification.

## **Lifetime battery (products using battery, only)**

KNOP calculates the life of the batteries. See specification in the user manual.

## **Intended use**

In general, the different variants of the medical device (transmitters/receiver system from KNOP Elektronik) are made as assistive aid for impaired/disabled patients to be able to call for assistance.

The various variants of the medical equipment consist of a combination of products (transmitter and receiver) designed to call for help to impaired/disabled patients; The transmitter-receiver system can be activated in different ways:

- For example, patients who actively do this and are aware that they are calling for help, such as e.g., people with walking difficulties who need help going to the toilet.
- Or disabled patients who are not aware that they are inadvertently putting pressure on the sendere.g. during a seizure.



- Or patients with intellectual and cognitive deficits who are not aware that a receiver is receiving a signal from their transmitter when they leave a house or room.

The system is not designed for critically or mentally ill persons.

### **General product description**

The products manufactured by KNOP elektronik consist of several variants of transmitters and receivers that can be combined with each other. In addition, these products are used in combination with positioning and repeater systems.

These products are medical devices intended to call for assistance and are used for disabled/mobile people, such as those with walking difficulties who need help going to the toilet. The transmitters are activated, for example, by a sound or by pressure (e.g. by pressure with a finger or a breath through the mouth). The recipient is supervised by health care personnel or lay people in private homes.

The system is not designed for critically or mentally ill people.

Part of product	Function in the product system
Transmitter	The transmitter can send the signal obtained from the patient to the receiver monitored by the healthcare personnel or lay person. The transmitter products can be activated by button, sound, blow or movements.
Receiver	The transmitters can be coded into all the receivers and to several receivers at the same time. Some receivers also have a summon button for calling assistance from their coworkers.
Repeatersystem	If it is needed to cover a more comprehensive and larger area a repeater system is used. The repeater system also gives an increased functionality as e.g. that an alarm automatically is received at first at the healthcare person closest to the client.
Positionsystem	If a sender is equipped with a position receiver it can be used in connection with a Position system. Not all product variants include position receivers. Typically, it is seen in connection with nursing homes and security for patients with dementia. The receivers can be portable or stationary.

### **Product description**

MR902 is the primary unit of our Master Repeat System.

Together with RP902, the MR902 will extend the range and functionality of the system and devices. Up to 14 repeaters can be added to a system.

As a stand-alone unit, it adds new functionality to the devices.

### **Installation**

MR902 only needs to be connected to a 230V socket. There are no cables between the units.



The product switches on automatically when it is connected to the mains.

#### **Light emitting diode - LED**

The LED is either green or red when switched on.

<b>RED:</b>	Error, not connected.
<b>RED flashing:</b>	Battery mode.
<b>GREEN:</b>	Connected, everything works fine.
<b>GREEN flashing:</b>	Confirmation.
<b>NO light:</b>	The device is switched off.



### Backup

MR902 has an internal backup battery to keep the system running in the event of a power failure. The battery lasts for a few hours. If necessary, the battery can be recharged (full charge takes about 24 hours). The battery is not intended to be used as a replacement. Contact your retailer for battery service.

### Commissioning

Before using MR902, the system must be configured for the current system.

**Important:** It is important that the clock is set on the MR902 as it is the system clock. Other products will be set after MR902 system clock. See main menu.

Note: The display language is in English only.

Setting up a system is quite simple.

- First determine the number of repeaters. In most cases, it is fine if you have no more than 25 meters between them. It depends on the building material, is it wood, concrete etc. You can check the connection level in the menu. (menu reports) We recommend above 3.0, otherwise you can just add an extra RP902.
- For example, if you have 5 repeaters, you have set the system size to 6, which is closest to the menu settings/system size.
- When you need to add repeaters to the system (menu repeaters / add new)
- Wait for about one minute per repeater. The system will automatically figure out how to make the connections between them. When all RP902 and MR902 units are ready, the green LED will light up. This is done every time the system is switched on. When the system is running normally, the system often checks the interconnection in the background.
- The basic system is now working. The installation of RP902 and MR902 is now complete.
- Receivers must now be added, this is done in the receiver menu when setting up the system or via KNOPtool when connecting to MR902. See the user manual for your receiver.
- KNOPtool can be used to configure the whole system as well as the receiver/transmitter. See the manual for installation of KNOPtool. This way you have a complete copy of the system setup.




## Start screen

The unit must be connected to a wall outlet before use.



The system clock normally appears on the display.

The system name "if coded" is shown at the top. By default, MR902.



*Set the clock !*


## Main menu


If you press  and  at the same time, either on the lock screen or on the home screen, it opens the main menu.


General information about the menu function:



## Operation

Press   to move the green arrow either up/down or right/left in the menu.

Press  to cancel the selection or exit the menu.

Press  to accept the current choice.

**Note:** After about 15 seconds of inactivity, the menu returns to the home screen.

## Main menu - Repeaters


In this menu it is possible to add, modify, delete and check the status of the repeater.



Select **Repeaters** and press  to confirm. Press  to exit.

**Add a new repeater, press "Add new".**

This means that a new repeater is added to the system.




exit.

Select **"Add new"** and press  to confirm. Press  to



Please wait, the MR902 will search for a new RP902 device.



To add the RP902 device found, press **Ent** to confirm. Press **Esc** to exit.


The display has now been added to RP902 unit. The LED on RP902 will also glow green to show that it is connected.



In this case, the added RP902 is shown on the top line as repeater "1".  
"1:5.0Y 9.9%", where 1=RP902 number, 5.0Y=antenna quality, 9.9%=rebroadcasts.

For a newly started system, the retransmissions will show a higher value due to the set transmissions at start-up, if needed you can select "Reset statistics in menu".

If you select a repeater by pressing **Ent** on the repeater you want to view, you will get a list with information about the repeater.



**ID:** The ID or serial number of the selected repeater.  
**Err:** Transmission error.  
**OK:** Call handled.  
**RX:** Receive level. Higher numbers are better, e.g. -30 dBm is better than -40 dBm.  
**Link:** Connected for example to a Master or another RP902.  
**Discon:** Number of interruptions.  
**Noise:** Noise level. Lower numbers are better, e.g. -112 dBm better than -100 dBm.



## Restored status.

This clears all status logs etc.



Select **"Reset status"** and press **Ent** to confirm.


## Edit name

Edit the name of the repeater. The name is not used in the system, but it is for your information. The name can be the location of the repeater.



Select **"Edit name"** and press **Ent**. Press **Esc** to exit.

Press **Esc** to return to the Repeater menu.



Press **↓** **↑** to change the letter.

Press **Ent** to select a letter.

Press **Esc** to choose to delete a letter.


Press **Esc** for about 3 seconds to close the menu.

To save the new name, select the symbol **♦** and press **Ent**.




## Delete

Remember to erase a repeater unit if it is not used or is replaced. An idle or replacement can slow down the system.




Select delete and press **Ent**. And press **Esc** to leave.



Pressing **Ent** confirms that it should be deleted. Press **Esc** to leave.

## Main menu - Receiver

In this menu you can delete a receiver from the system and check its status.



Select the desired recipient and press **Ent** to confirm. Press **Esc** to exit

<b>Product:</b>	Product type e.g. RX901B
<b>ID:</b>	Product ID/serial number, e.g. 63123.
<b>Name:</b>	Product name e.g. RX901B
<b>SW ver.:</b>	Software version of the product, e.g. 3.4.17.
<b>Online:</b>	Status, online or offline.
<b>MR902:</b>	MR902 reception status in the product e.g. -33dBm.
<b>RP901-1:</b>	Reception status of RP902 In the product e.g. -53dBm.
<b>Batt:</b>	Battery level, either voltage or percentage, e.g. 20%.

**Tip:** If you name your recipients, it's easy to find them in a list like this, otherwise it's a list with the same name.

## Delete

Remember to remove a receiver device if it is not in use or replace it if it is inactive or offline.



Select **"Delete"** and press **Ent**. Press **Esc** to leave.



Press **Ent** to confirm the deletion... Press **Esc** to exit.



## Main menu - Settings

This menu contains general information about system and clock settings, etc.

Main menu	Settings	Settings	Settings
Repeaters Receivers → Settings	ID: 58961 Ver. 1.6.4 Time: 15 System: 6 ↓ Set clock	↑ System: 6 → Set clock Set reply time Edit name ↓ Edit system	↑ Edit name Edit system System test Reset all → Power off

Select "**Settings**" and press **Ent** to confirm. Press **Esc** to exit.

**ID:** ID is the normal serial number of an MR902.  
**SW ver.:** Software version (firmware) of an MR902.  
**Reply time:** The reply time can also be set.  
**System size:** Size of the system, maximum number of repeaters.

## Set the clock

It is important that the clock in MR902 is set correctly. If you use timers in another product that is connected or used in the system, these timers will work according to this clock. The clock in other products is usually also synchronized with the system clock.

Main menu	Settings
Repeaters Receivers → Settings	↑ System: 6 → Set clock Set reply time Edit name ↓ Edit system

Select "**Set clock**" and press **Ent** to confirm. Press **Esc** to exit.

Set clock	Set clock	Set clock	Set clock	Set clock
2021-02-17 08:56 Year: 2021	2021-02-17 08:56 Month: 02	2021-02-17 08:56 Hour: 08	2021-02-17 08:56 Day: 17	2021-02-17 08:56 Minute: 56

Press **↓** **↑** to change the setting.


Press **Ent** to accept and move on.

Press **Esc** to leave.



## Set response time

This indicates the time between alarms on the receivers. When the system receives an alarm, MR902 estimates the location of the nearest receiver and then calls that unit. If the call is not answered or overheard within this time, the call goes to the next receiver. Normally, 15 seconds is a typical time to accept the call. You can change the time from 5 to 30 seconds.



Select **"Reply time"** and press **Ent** to confirm. Press **Esc** to leave.




Press **↓** **↑** to change the setting.

Press **Ent** to approve.


Press **Esc** to leave.

## Edit name

Change the system name. The name is displayed on the main screen.



Select **"Edit name"** and press **Ent**. Press **Esc** to exit.



Press to change the letter.

Press to select a letter.


Press to choose to delete the last letter.

Press to exit the menu. It takes 3 seconds.


To save the new name, select the and press .

## Size of the system

Set the size of the system as close as the options allow, but the number should be greater than or equal to the number of repeaters. Response time will be faster the lower the value you choose.



Select "**System size**" and press . Press to exit.



Press to resize the system.

Press to accept the new choice.

Press to leave.



## Systemtest

System testing is a troubleshooting tool to update statistics as quickly as possible. This function should normally be turned off. If you turn it on, it is important to turn it off again.

Main menu	Settings	System test
Repeaters Receivers → Settings	↑ Set reply time Edit name Edit system → System test ↓ Reset all	→ Reset stat. RP test: Off Rx test: Off

Select **Reset stat.** and press **Ent** to reset all statistics.

Select **RP test** and press **Ent** to activate or deactivate the repeater test.

Select **RX test** and press **Ent** to activate or deactivate the repeater test.

Press **Esc** to exit.

## Reset all

**CAUTION!! This will delete everything and be removed from all devices.**

Main menu	Settings
Repeaters Receivers → Settings	↑ Edit name Edit system System test → Reset all Power off

Select **"Reset all"** and press **Ent**. Press **Esc** to exit.

Select **Ent** to accept to **RESET ALL!**



Select **Esc** to leave.






## Turn off

This shuts down MR902 and the system. All calls are no longer handled by the system.

**IMPORTANT:** Before selecting "**Power off**", MR902 must be unplugged from the power outlet.



Select "**Power off**" and press **Ent** and press **Esc** to exit.



Select **Ent** to accept to **DELETE EVERYTHING!**

Select **Esc** to leave.



## ***Maintenance***

### **Updating software**

It is possible to update the software of MR902. For a more detailed description please check the KNOP tool.

### ***Cleaning***

The product can be cleaned with a moist cloth or a dry cloth.

### ***Spare parts and accessories***

Spare parts and accessories can be ordered from Knop Elektronik A/S at [www.knop.dk](http://www.knop.dk).



## Technical data

<b>Frequency HF:</b>	869.2125 MHz
<b>Range:</b>	Up to 1500 m between MR902/RP902 and RP in free sight. <sup>(1)</sup>
<b>Range in building:</b>	Recommended distance between repeater units <25m
<b>Power supply:</b>	230VAC
<b>Connection plug:</b>	2-pole EURO
<b>Battery:</b>	Lithium battery, LP853048 3.7V 1200mA
<b>Battery time:</b>	Approximately 4 hours
<b>Charge:</b>	Approximately 1 day
<b>Surrounding environment:</b>	Indoor use. ≤ 90% non-condensing
<b>Ambient temperature:</b>	0°C to +40°C
<b>Cabinet type:</b>	Color RAL9001 ABS
<b>Cabinet dimensions:</b>	Ø110mmx52mm
<b>Desnity:</b>	IP40
<b>Weight:</b>	135g















The right to make changes is reserved

All rights reserved.

© **KNOP** Elektronik A/S

<sup>(1)</sup> Measured outdoors with an unobstructed view between transmitter and receiver. In buildings, the range is reduced.

## Explanation of symbols and approvals

	<b>This product complies with:</b> Directive 2017/745/EU      MDR Directive 1907/2006/EU      REACH Directive 2011/65/EU      RoHS Directive 2012/19/EU      WEEE ISO 14971:2019      Risk Management for Medical devices EN 301 498-1 V2.2.3      Electro Magnetic Compatibility EN 301 489-3 V3.1.1      Electro Magnetic Compatibility EN 50130-4:2011 + 2014      Immunity alarm systems EN 300 220-1 V3.1.1      Short Range Devices EN 300 220-2 V3.1.1      Short Range Devices EN 300 220-3 V2.1.1      Short Range Devices EN 60601-1-2:2014 + 2015      Electromagnetic compatibility (Medical) EN 62368-1:2020      Electrical safety EN 50581:2012      Hazardous substances
	The product must not be used if the packaging is damaged.
	<b>Medical Device</b> Class 1, rule 1
	<b>Manufacturer</b> KNOP Elektronik A/S, Fabriksvej 20, DK-7600 Struer, Denmark
	Read the manual(s) before installation and commissioning at <a href="http://www.knop.dk">www.knop.dk</a> .
	Interference may occur in an environment with equipment marked with this symbol.
	Must be protected against liquids.
	0 °C to +40 °C, temperature limit for transport/storage and use.
	The product must not be disposed of with normal household waste.
	<b>Single Registration Number</b> DK-MF-000025631
	<b>Unique Device Identifier</b> 05744002850335
	<b>Product reference</b> MR902 Master Repeater
	<b>Serial number</b> Placed on the product
	<b>Ingress Protection Code</b> IP40