

User manual

RX901B





Table of contents:

Warning and safety instructions:	3
Product description	5
Commissioning	6
Reception of alarms	8
Emergency call	9
Batteryalarm	10
Clients and transmitters	11
Main menu	12
Settings and encoding:	14
Maintenance	23
Cleaning	23
Spare parts and accessories	23
Check the coverage area	24
Technical data	25
Explanation of symbols and approvals	26

Warning and safety instructions



- The battery in this product complies with the relevant IEC safety standards for batteries.
- Please read the safety instructions on the power supply before use.
- The power supply unit must meet the relevant IEC safety standards.
- When replacing the power supply unit or the charging plate, original and/or approved units must be used.
- The battery must not be replaced, and the product must be sent for repair.
- Please read the intended use before putting it into service.
- The product is not water-resistant unless noted in product specification.

Contraindications

In general, the product cannot be used if:

- The disabled client/patient is mentally ill.
- The disabled client/patient is critically ill.
- The client/patient is unable to use the transmitters.

Lifetime after mounting of the device

The battery, if applies, must be replaced when the battery replacement information appears on the display and according to the user manual.

Lifetime is evaluated in relation to the pressure button. This is made according to the supplier's specification.

Lifetime battery (products using battery, only)

KNOP calculates the life of the batteries. See specification in the user manual.

Intended use

In general, the different variants of the medical device (transmitters/receiver system from KNOP Elektronik) are made as assistive aid for impaired/disabled patients to be able to call for assistance.

The various variants of the medical equipment consist of a combination of products (transmitter and receiver) designed to call for help to impaired/disabled patients; The transmitter-receiver system can be activated in different ways:

- For example, patients who actively do this and are aware that they are calling for help, such as e.g., people with walking difficulties who need help going to the toilet.



- Or disabled patients who are not aware that they are inadvertently putting pressure on the sender e.g. during a seizure.
- Or patients with intellectual and cognitive deficits who are not aware that a receiver is receiving a signal from their transmitter when they leave a house or room.

The system is not designed for critically or mentally ill persons.

General product description

The products manufactured by KNOP elektronik consist of several variants of transmitters and receivers that can be combined with each other. In addition, these products are used in combination with positioning and repeater systems.

These products are medical devices intended to call for assistance and are used for disabled/mobile people, such as those with walking difficulties who need help going to the toilet. The transmitters are activated, for example, by a sound or by pressure (e.g. by pressure with a finger or a breath through the mouth). The recipient is supervised by health care personnel or lay people in private homes.

The system is not designed for critically or mentally ill people.

Part of product	Function in the product system
Transmitter	The transmitter can send the signal obtained from the patient to the receiver monitored by the healthcare personnel or lay person. The transmitter products can be activated by button, sound, blow or movements.
Receiver	The transmitters can be coded into all the receivers and to several receivers at the same time. Some receivers also have a summon button for calling assistance from their coworkers.
Repeatersystem	If it is needed to cover a more comprehensive and larger area a repeater system is used. The repeater system also gives an increased functionality as e.g. that an alarm automatically is received at first at the healthcare person closest to the client.
Positionsystem	If a sender is equipped with a position receiver it can be used in connection with a Position system. Not all product variants include position receivers. Typically, it is seen in connection with nursing homes and security for patients with dementia. The receivers can be portable or stationary.



Product description

RX901B is a receiver designed to work with KNOP 90X series wireless devices.

The unit can be worn in a pocket or on a necklace around the neck.

Up to 250 transmitters can be encoded and assigned a name of your choice.

The last 50 calls are stored in a list where you can see the time of the call, where the call came from, and whether it was acknowledged.

With an RX901B, it is possible, for example, to divide patients (residents) in a care institution into groups, e.g. in connection with shift schedules, where the night shift has to serve more patients than the day shift.

RX901B is equipped with a position reader, so it can be used in the KNOP positioning system.

With RX901B it is possible to make an emergency call to another colleague with a single touch.

Operation is menu-structured in a graphic display and controlled by the keyboard.

RX901B is equipped with a rechargeable battery and wireless charging.

The product is waterproof (IP65).

Commissioning

It is recommended that RX901B be charged to 100% before use.

The supplied wireless charger connects to the adapter via the USB cable.

RX901B is placed on top of the charging plate and the RX901B's display lights up, showing the battery's charge level in percent.

If the device is fully discharged, it may take up to 5 minutes for the display to light up.

Package contents:





RX901B in the charger:



When the product is switched on for the first time, the desired language must be selected.



Press   to move the green arrow to the desired language.

Press  to select a language.



Home screen

The first time the device is to be switched on after unpacking, it should be placed on the charger.

When RX901B is not connected to the charger, the screen will turn off approximately 10 seconds after the last operation.

Pressing any key turns on the display and shows the time and date.

The user's name is displayed at the top, otherwise, it is coded "RX901B".



Home screen options:



The name of e.g. a product or a location can be changed to the desired name.

The current time and date are displayed in white text.

The clock must be set.

Battery mode:



Battery OK.



Charging recommended.



The battery is almost empty and needs charging.

Alarm sound:



The alarm sound is on.



The alarm sound is switched off.

Vibrator:



The vibrator is activated.



The vibrator is switched off.

System - displayed only if the system is connected:



The system is connected.



Connecting to the system.



Interruption of the system signal.



The system signal has been lost.

Press  to see the call list.



Reception of alarms

When the receiver receives an alarm from a coded transmitter, the name, type and time of the reception and any position are displayed.



Client's name
Optional text if encoded.
Type of transmitter.
Positioning for the alarm.
Button/function type.

The individual texts are only displayed if they are encoded or known.


Press  to confirm the alarm.

Press  to reject.

Note: In the simple system, the alarm call is handled in the same way, regardless of whether it is accepted or rejected, instead of the KNOP repeater system (MR90X/RP90X) where the alarm is automatically forwarded to other receivers.

If the alarm is not answered, short sounds are heard approximately every five minutes.

When a message is received from a transmitter, e.g. with a low battery, the message below is displayed.



Press  to confirm the alarm.

Note: "Last call" displays the "low battery" message for the transmitter where the battery needs to be replaced.



Emergency call

It is possible to send an emergency call (or alarm) from the recipient to e.g. other users to indicate "I need help" at a patient or as an assault alarm.

The emergency call is activated by pressing the KNOP logo.

When the emergency call is sent, the receiver vibrates once and when the emergency call is received by another user, the receiver vibrates three times.




Emergency call	Emergency call	Emergency call
Transmit Emergency call!	Transmit Emergency call!	Transmit Emergency call!
	OK	Error!
Transmitting	Accepted	Failed



Batteryalarm

The receiver indicates the battery status on the display and in case of discharged battery with sound indications.

Battery mode:

-  Battery OK.
-  Charging recommended.
-  The battery is almost empty and needs charging.

If the symbol changes to red, put the device in the charger immediately.

At the same time, a few short beeps will be heard about every five minutes, indicating low battery, while the screen below is displayed:





Clients and transmitters

Before you can encode a transmitter in the receiver, you must create a client to which the transmitter belongs.

A client can be a person, a room, or something else. For example, the client can be **Sven** (used as an example).

Our examples our client is called **Sven** can have one or more transmitters associated with him.

For example, **Sven** (our client) has a wristband transmitter, a bed alarm, a motion sensor, and other products.

When **Sven** activates an alarm from his products, the display shows that the alarm comes from **Sven** and which product is the sender, so that the cause of the alarm is known from the beginning.

This of course assumes that names etc. are encoded in the desired way.

Clients created under Settings / Client / New client:



Main menu	Settings	Clients
↑ Hibernate Select Group Vibrator: On Sound: Low → Settings	Ver: 3.8.1 ID: 123456 Clients → Edit group ↓ Set clock	→ Add client Sven All TX's.

Add under Setting / Client / Client (Sven) / Add TX unit:


Main menu	Settings	Clients	Sven
↑ Hibernate Select Group Vibrator: On Sound: Low → Settings	Ver: 3.8.1 ID: 123456 Clients → Edit group ↓ Set clock	Add client → Sven All TX's.	↑ Add TX unit Select group Client name → Client text Delete



Main menu

If you press  and  at the same time, either when the screen is off or on the home screen, the main menu opens.

The green arrow on the left shows which line in the menu can be selected.



Operation:

Press  to move the green arrow either up or to the right in the menu.


Press  to move the green arrow either down or left in the menu.

Press  to leave the menu.

Press  to select the menu item with the green arrow next to it.

NOTE: After about 60 seconds of inactivity, the menu closes automatically.


Latest alarms:




  Scroll between alarms.

GREEN Means that the alarm was accepted.
RED Means that the alarm was rejected.
YELLOW Means the alarm went unanswered.

See more information about the alarms:



 See the client's position and which transmitter the alarm is coming from.

  Scroll through multiple alarms.




When an alarm is received from a transmitter that has, for example, a low battery, the message is displayed.

  to scroll through other alarm calls.



Delete the call list:



Press **Ent** to confirm.


Suspend/Hibernate receiver:



Ent Confirm or **Esc** to reject.

The receiver is switched on by pressing **Ent**.

Group selection:



Ent Select or deselect a group.


Esc Return to menu.

↓ **↑** Moves arrow to the group you want to change.

1 Group selected

1 Group deselected


Vibrator function:



Ent Activate or deactivate the vibrator function.

Esc Return to menu.

Sound setting:



Ent Choose between **High**, **Low** or **Off** sound when receiving alarm calls.



Settings and encoding:

Main menu	Settings	Settings	Settings
↑ Hibernate Select Group Vibrator: Off Sound: High → Settings	Ver: 3.8.1 ID: 064067 → Clients Edit group ↓ Set clock	↑ Set clock Click: Low Product name → Passcode ↓ Outgoing	↑ Outgoing System menu Positions Language → Turn off

Ver: 3.3.14

The software version of the product.

ID: 58651

Serial number of the product.

Add new client:

Settings	Clients
Ver: 3.8.1 ID: 123456 → Clients Edit group ↓ Set clock	→ Add client Sven All TX's.

Select **Add client** or select client and press **Ent**.

It will ask for the client's name and optional text.
Remember to add TX devices.

Change client:

Sven	Clients	Sven
↑ Add TX unit Select group Client name → Client text Delete	Add client → Sven All TX's.	Room 10 → TX901 Add TX unit Select group Client name ↓

Sven
Room 10
→ TX901


Client name.

Optional text, e.g., name of room.


Product already encoded on the client.



Add transmitters:



Now activate the desired transmitter.




Accept new transmitter.



If the transmitter is already encoded, the menu below appears, allowing you to move a transmitter between clients.


Group selection of transmitter:



- Select or deselect a group.
- Return to menu.
- Moves arrow to the group you want to change.

- Group selected
- Group deselected

Change the client's name:



- Choose letter.
- Select letter.
- Delete last letter.

To save the text select symbol and press .

Add/change client optional text:







- Choose letter.
- Select letter.
- Delete last letter.

To save the text select symbol and press .



Delete client:

Sven	Sven
↑ Add TX unit	Delete:
Select group	Sven
Client name	Confirm?
Client text	 
→ Delete	

 Confirm delete.
 Return to menu.



Change encoded transmitter:

Settings	Clients	Sven	Sven	Sven
Ver: 3.8.0 ID: 48 → Clients Group ↓ Set clock	Add client → Sven Allan All Tx's.	Room 10 → TX901 Add TX unit Select group ↓ Client name	↑ Optional text Alarm signal Call texts → Call: Normal ↓ Timer setup	↑ Alarm signal Call texts Call: Normal Timer setup → Delete

TX901
ID: 60372
Ver.: 1.0
Batt.: 40%
Signal 1

Product name of the transmitter.
Transmitter ID.
The software version of the transmitter.
Battery level in the transmitter (if applies).
Selected audio signal when alarm is received.

Change name of the transmitter:

Sven	Optional text
↑ Batt.: 40% Signal 1 → Optional text Alarm signal ↓ Call texts	Save? ◆ 456789 a b c d Enter TX optional text

↓ ↑ Choose letter.

Ent Select letter.

Esc Delete last letter.

To save the text select symbol ◆ and press Ent.

Change alarm signal:

Sven	Sven
↑ Batt.: 40% Signal 1 Optional text → Alarm signal ↓ Call texts	TX901 Signal 1 ↓ 1 2 3 4 5

↓ ↑ Scroll between audio signals.

Ent Select audio signal.

Change call text:

Sven	Edit text	1. Green
↑ Signal 1 Optional text Alarm signal → Call texts ↓ Call: Normal	→ 1. Green 2. Red 3. Yellow 4. Blue ↓ KNOP button	Alarm button → 1. Green 2. Red 3. Yellow ↓ 4. Blue

This is a TX901 example.
It is possible to change the shown text
for a button etc.
In KNOPtool used defined texts can be
made.



Reset call text:

Sven	Edit text	Call texts
↑ Signal 1	↑ Extern 2	Reset texts
Optional text	Extern 3	
Alarm signal	Extern 4	Confirm?
→ Call texts	Low battery	 
↓ Call: Normal	→ Reset texts	

Restore call texts to factory setup.

Accepted alarm types:



Sven	Sven
↑ Optional text	↑ Optional text
Alarm signal	Alarm signal
Call texts	Call texts
→ Call: Normal	→ Call: All
↓ Timer setup	↓ Timer setup

 Toggle between **Normal** and **All**.



All accepts both types of alarms Normal and Emergency calls.

Transmitter activation timer:

Sven	Sven
↑ Alarm signal	→ Timer: Off
Call texts	Start 1 00:00
Call: All	Stop 1 00:00
→ Timer setup	Start 2 00:00
Delete	Stop 2 00:00

 Toggle between **On** and **Off**.  Change hour/minutes. Change selection.

Remove unit from client:

Sven	Sven
↑ Alarm signal	Delete:
Call texts	TX901
Call: All	Confirm?
Timer setup	 
→ Delete	

 Delete confirm. Return to menu.



Group settings:

Settings	Setup group
Ver: 3.8.1	Group: 1
ID: 123456	
Clients	
→ Edit group	↓
↓ Set clock	1 2 3 4 5 6 7 8

- ↓ ↑ Move between groups.
- Ent Select group.

Group activation timer on/off:

Group: 1	Group: 1
→ Timer: Off	→ Timer: On
Group name	Group name
Timer setup	Timer setup

- Ent Toggle between **On** and **Off**.
- Esc Return to menu.

Don't forget to set the time in "Settings - Set clock".

Change group name:

Group: 1	Group: 1
Timer: On	Save?
→ Group name	◆ 4 5 6 7 8 9 a b c d
Timer setup	Gruppe: 1
	Enter
	group name

- ↓ ↑ Choose letter.
- Ent Select letter.
- Esc Delete last letter.

To save the text select symbol ◆ and press Ent.

Change the timer:

Group: 1	Group: 1
Timer: On	Start 07:00
Group name	Stop 15:00
→ Timer setup	

- ↓ ↑ Change hours or minutes.
- Ent Select items.
- Esc Return to menu.

Set clock:

Settings	Set clock
↑ ID: 123456	2022-10-27
Clients	11:21
Edit group	Automatic
→ Set clock	Year: 2022
↓ Click: Low	

- ↓ ↑ Change item.
- Ent Scroll between **Year**, **Months**, **Time** and **Date**.

Note: If the receiver receives a clock update from a KNOP product or USB901, the clock is automatically set accordingly.



Button click sound:

Settings	Settings
↑ Clients	↑ Clients
Edit group	Edit group
Set clock	Set clock
→ Click: Low	→ Click: High
↓ Product name	↓ Product name

Ent Toggle between **Low**, **High** and **Off**.

Esc Return to menu

Change the product (user) name:

Settings	Name
↑ Edit group	Save?
Set clock	◆ 456789 a b c d
Click: High	RX901
→ Product name	Enter
↓ Passcode	product name

↓ **↑** Choose letter.

Ent Select letter.

Esc Delete last letter.

To save the text select symbol ◆ and press **Ent**.

Change passcode:

Settings	Passcode
↑ Set clock	New code
Click: Low	
Product name	
→ Passcode	Code: 1234
↓ Outgoing	

↑ **↓** Select digits (0 to 9).

Esc Delete last digit.

The code must always be 4 digits.

Function of the KNOP button:

Settings	Outgoing
↑ Click: Low	→ Call: All
Product name	
Passcode	
→ Outgoing	
↓ System menu	

Ent Change function, choose between **All**, **Normal** or **Off**.


Connecting and changing system code:

Settings	System menu
↑ Product name	→ System: Off
Passcode	ID: 0
Outgoing	
→ System menu	
↓ Positions	

Ent Select either system **Off** or **On**.



Set the system ID:




Press to change the value (0 to 9).

Press to select a number.


Press to delete the last digit.

Positions:




Current position.
Already encoded.

View position:



Shows the current position of where you are.


Add position:



Save new position
 Return to menu.

Note that the receiver must be within range of a position field.

Change position name:



Choose letter.

Select letter.

Delete last letter.

To save the text select symbol ♦ and press .



Remove position:

Edit position	Edit position
Pos.: 255	Pos.: 255
Door	Door
Name	Delete?
→ Delete	<input type="checkbox"/> <input type="checkbox"/>

- ☐ Delete position.
☐ Return to menu.

Languages:

Settings	Language
↑ Outgoing	→ Dansk
System menu	Svenska
Positions	Norsk
→ Language	Suomalainen
Turn off	English

- ☐ ☐ choose language.
☐ Accept the chosen language.
☐ Return to menu.

Turn off the product:

Settings	Turn off
↑ Outgoing	
System menu	
Positions	
Language	Confirm?
→ Turn off	<input type="checkbox"/> <input type="checkbox"/>

- ☐ Confirm **Turn off**.
☐ Return to menu.

Note: The clock and date are reset.

To turn on the receiver, place the receiver on the charger.



Maintenance

Software update of RX901B:

It is possible to update the software in RX901B, as well as to download and transfer settings with the KNOPtool. For further descriptions please check the KNOPtool manual.

Cleaning

The product can be cleaned with a moist cloth or wipe.

Spare parts and accessories

Spare parts and accessories can be ordered from Knop Elektronik A/S at www.knop.dk.



Check the coverage area

One person activates a transmitter (e.g. TX901) at short intervals, while another person systematically walks around the area and marks on a sketch the building and the area where there is coverage. The sketch is placed in the office so that everyone who has to manage the system can see where there is coverage. Outside the coverage area, alarms cannot be received.

The coverage area can be extended using the KNOP repeater system.

Contact your retailer for more information.



Technical data

Battery:	1 pc. Li-Ion 3.7V 450 mAh
Operating time with battery:	Approximately 3 days with 50 alarm calls per day.
Charging:	Wireless, approx. 4 hours.
Power supply unit:	230VAC adapter USB 5VDC/1.0A
Frequency HF:	869.2125 MHz
Range:	Up to 1500 m free range for an RX901B ⁽¹⁾
Display:	Graphical LCD
Menu texts:	5 languages: Danish, Finnish, English, Norwegian, and Swedish.
Log function:	Saves the last 50 alarm calls.
Memory:	250 transmitters 150 positions POS901
Surrounding environment:	Indoor use. < 90% humidity without condensation.
Ambient temperature:	0 - 40°C
Cabinet type:	Color RAL9003 ABS
Cabinet dimensions:	45 x 90 x 16 mm. Incl. neck strap with safety lock that breaks when you pull it.
Density:	IP65
Weight incl. batteries:	56g















The right to make changes is reserved.

All rights reserved.

© **KNOP** Elektronik A/S

⁽¹⁾ Measured outdoors with an unobstructed view between transmitter and receiver. In buildings, the range is reduced.

Explanation of symbols and approvals

	This product complies with: Directive 2017/745/EU MDR Directive 1907/2006/EU REACH Directive 2011/65/EU RoHS Directive 2012/19/EU WEEE ISO 14971:2019 Risk Management for Medical devices EN 301 498-1 V2.2.3 Electro Magnetic Compatibility EN 301 489-3 V3.1.1 Electro Magnetic Compatibility EN 50130-4:2011 + 2014 Immunity alarm systems EN 300 220-1 V3.1.1 Short Range Devices EN 300 220-2 V3.1.1 Short Range Devices EN 300 220-3 V2.1.1 Short Range Devices EN 60601-1-2:2014 + 2015 Electromagnetic compatibility (Medical) EN 62368-1:2020 Electrical safety EN 50581:2012 Hazardous substances
	The product must not be used if the packaging is damaged.
	Medical Device Class 1, rule 1
	Manufacturer KNOP Elektronik A/S, Fabriksvej 20, DK-7600 Struer, Denmark
	Read the manual(s) before installation and commissioning at www.knop.dk .
	Interference may occur in an environment with equipment marked with this symbol.
	Must be protected against liquids.
	0 °C to +40 °C, temperature limit for transport/storage and use.
	The product must not be disposed of with normal household waste.
	Single Registration Number DK-MF-000025631
	Unique Device Identifier 05744002850298
	Product reference RX901B Receiver
	Serial Number Placed on the product
	Ingress Protection kode IP65