

User manual

KK100





Table of contents

Warning and safety regulations.....	3
Intended use.....	4
Commissioning.....	4
Connection.....	4
Service.....	5
Location.....	5
Cleaning.....	5
Spare parts and accessories.....	5
Technical data.....	6
Explanation of symbols and approvals.....	7



Warning and safety regulations



- Read the intended use before putting into service.



Intended use

The product is designed to be placed close to the user's head, for example, so that the alarm can be activated by turning the user's head towards KK100.

KK100 can be used when the user cannot activate an alarm with a pull cord or push button.

KK100 must be used in conjunction with a KNOP Elektronik transmitter or existing call systems.

Commissioning

After installation, it is important that the user and staff test the product.

Staff should receive an alarm on their alarm device/receiver.


Please refer to the user manual of the transmitter or paging system used.

Connection

KK100 has a 3m cord with 3.5mm jack.

Potential-free contact between frame and tip.

Max. 36V/100mA.





Service

When the user activates KK100, an alarm is sent via the connected transmitter or pager.

Location

KK100 shall be positioned so that the user can easily activate the ball on the switch. The cheek switch can be mounted on a gooseneck, for example. Please contact your dealer for further information.

Cleaning

The product can be cleaned with a damp cloth or a disinfected wet wipe.

Spare parts and accessories

Spare parts and accessories can be ordered on our website www.knop.dk



Technical data

Surrounding environment:	Indoor use \leq 90% non-condensing
Ambient temperature:	0°C to +40°C
Enclosure type:	White polyoxymethylene
Enclosure dimensions:	Ø: 43mm L: 106mm
Thread for mounting:	3/8"
Density:	IP40
Weight:	90g















The right to make changes is reserved.

All rights reserved.

KNOP Elektronik A/S



Explanation of symbols and approvals

	<p>This product complies with these directives and standards:</p> <table border="0"> <tr> <td>Directive 2017/745/EU</td> <td>MDR</td> </tr> <tr> <td>Directive 1907/2006/EU</td> <td>REACH</td> </tr> <tr> <td>Directive 2011/65/EU</td> <td>RoHS</td> </tr> <tr> <td>ISO 14971:2019</td> <td>Risk Management for Medical devices</td> </tr> <tr> <td>EN 301 498-1 V2.2.3</td> <td>Electro Magnetic Compatibility</td> </tr> <tr> <td>EN 301 489-3 V3.1.1</td> <td>Electro Magnetic Compatibility</td> </tr> <tr> <td>EN 50130-4:2011 + 2014</td> <td>Immunity alarm systems</td> </tr> <tr> <td>EN 300 220-1 V3.1.1</td> <td>Short Range Devices</td> </tr> <tr> <td>EN 300 220-2 V3.1.1</td> <td>Short Range Devices</td> </tr> <tr> <td>EN 300 220-3 V2.1.1</td> <td>Short Range Devices</td> </tr> <tr> <td>EN 60601-1-2:2014 + 2015</td> <td>Electromagnetic compatibility (Medical)</td> </tr> <tr> <td>EN 62368-1:2020</td> <td>Electrical safety</td> </tr> <tr> <td>EN 50581:2012</td> <td>Hazardous substances</td> </tr> </table>	Directive 2017/745/EU	MDR	Directive 1907/2006/EU	REACH	Directive 2011/65/EU	RoHS	ISO 14971:2019	Risk Management for Medical devices	EN 301 498-1 V2.2.3	Electro Magnetic Compatibility	EN 301 489-3 V3.1.1	Electro Magnetic Compatibility	EN 50130-4:2011 + 2014	Immunity alarm systems	EN 300 220-1 V3.1.1	Short Range Devices	EN 300 220-2 V3.1.1	Short Range Devices	EN 300 220-3 V2.1.1	Short Range Devices	EN 60601-1-2:2014 + 2015	Electromagnetic compatibility (Medical)	EN 62368-1:2020	Electrical safety	EN 50581:2012	Hazardous substances
Directive 2017/745/EU	MDR																										
Directive 1907/2006/EU	REACH																										
Directive 2011/65/EU	RoHS																										
ISO 14971:2019	Risk Management for Medical devices																										
EN 301 498-1 V2.2.3	Electro Magnetic Compatibility																										
EN 301 489-3 V3.1.1	Electro Magnetic Compatibility																										
EN 50130-4:2011 + 2014	Immunity alarm systems																										
EN 300 220-1 V3.1.1	Short Range Devices																										
EN 300 220-2 V3.1.1	Short Range Devices																										
EN 300 220-3 V2.1.1	Short Range Devices																										
EN 60601-1-2:2014 + 2015	Electromagnetic compatibility (Medical)																										
EN 62368-1:2020	Electrical safety																										
EN 50581:2012	Hazardous substances																										
	<p>The product shall not be used if the packaging is damaged.</p>																										
	<p>Medical Device Class 1</p>																										
	<p>Manufacturer KNOP Elektronik A/S, Fabriksvej 20, DK-7600 Struer</p>																										
	<p>Read the manual(s) before installation and commissioning.</p>																										
	<p>Interference may occur in an environment with equipment marked with this symbol.</p>																										
	<p>Must be protected against liquids.</p>																										
	<p>0 °C to +40 °C, temperature limit for transport/storage and use.</p>																										
	<p>The product must not be disposed of with normal household waste.</p>																										
	<p>Single Registration Number DK-MF-000025631</p>																										
	<p>Unique Device Identifier KK100: 05744002852049</p>																										
	<p>Product reference/product number/warning text. KK100: Cheek switch</p>																										
	<p>Serial number Placed on the product</p>																										
	<p>Ingress Protection code IP40</p>																										