

# User manual

# RX901C





## ***Table of contents***

Warnings and safety instructions	3
Product description	5
Connections	6
Commissioning	7
Start screen	8
Reception of alarms	9
Batteryalarm	11
Clients and transmitters	12
Main menu	13
Settings and encoding:	15
Maintenance	24
Cleaning	24
Spare parts and accessories	24
Technical data	25
Explanation of symbols and approvals	26

## **Warnings and safety instructions**



- The battery in this product complies with the relevant IEC safety standards for batteries.
- Do not consume or eat the battery, risk of chemical burns.
- Keep new and used batteries away from children.
- If the battery compartment cannot be closed properly, discard the product and keep it away from children.
- Read the intended use before use.
- Batteries should be checked regularly.
- The product is not water-resistant unless noted in product specification.

### **Contraindications**

In general, the product cannot be used if:

- The disabled client/patient is mentally ill.
- The disabled client/patient is critically ill.
- The client/patient is unable to use the transmitters.

### **Lifetime after mounting of the device**

The battery, if applies, must be replaced when the battery replacement information appears on the display and according to the user manual.

Lifetime is evaluated in relation to the pressure button. This is made according to the supplier's specification.

### **Lifetime battery (products using battery, only)**

KNOP calculates the life of the batteries. See specification in the user manual.

### **Intended use**

In general, the different variants of the medical device (transmitters/receiver system from KNOP Elektronik) are made as assistive aid for impaired/disabled patients to be able to call for assistance.

The various variants of the medical equipment consist of a combination of products (transmitter and receiver) designed to call for help to impaired/disabled patients; The transmitter-receiver system can be activated in different ways:

- For example, patients who actively do this and are aware that they are calling for help, such as e.g., people with walking difficulties who need help going to the toilet.



- Or disabled patients who are not aware that they are inadvertently putting pressure on the sender e.g. during a seizure.
- Or patients with intellectual and cognitive deficits who are not aware that a receiver is receiving a signal from their transmitter when they leave a house or room.

The system is not designed for critically or mentally ill persons.

### **General product description**

The products manufactured by KNOP elektronik consist of several variants of transmitters and receivers that can be combined with each other. In addition, these products are used in combination with positioning and repeater systems.

These products are medical devices intended to call for assistance and are used for disabled/mobile people, such as those with walking difficulties who need help going to the toilet. The transmitters are activated, for example, by a sound or by pressure (e.g. by pressure with a finger or a breath through the mouth). The recipient is supervised by health care personnel or lay people in private homes.

The system is not designed for critically or mentally ill people.

<b>Part of product</b>	<b>Function in the product system</b>
Transmitter	The transmitter can send the signal obtained from the patient to the receiver monitored by the healthcare personnel or lay person. The transmitter products can be activated by button, sound, blow or movements.
Receiver	The transmitters can be coded into all the receivers and to several receivers at the same time. Some receivers also have a summon button for calling assistance from their coworkers.
Repeatersystem	If it is needed to cover a more comprehensive and larger area a repeater system is used. The repeater system also gives an increased functionality as e.g. that an alarm automatically is received at first at the healthcare person closest to the client.
Positionsystem	If a sender is equipped with a position receiver it can be used in connection with a Position system. Not all product variants include position receivers. Typically, it is seen in connection with nursing homes and security for patients with dementia. The receivers can be portable or stationary.



## **Product description**

RX901C is a receiver in the KNOP 901 series.

250 transmitters can be encoded and assigned an optional name.

The last 50 calls are stored in a list where you can see the time of the call, where the call came from, and whether it has been acknowledged.

With RX901C it is possible to divide the residents of a care centre into groups, e.g. in connection with shift schedules, where the night shift has to serve more residents than the day shift.

RX901C is equipped with a position reader, so it can be used in KNOP's positioning system.

With RX901C it is possible to make an emergency call to another colleague with a single touch.

RX901C has several outputs that can be connected to external equipment, such as another pager, a phone transmitter with a call to the mobile phone, or a bed vibrator that can be placed under the pillow when you go to bed.

The controls are menu-driven with a graphical LCD screen and a click pad.


RX901C is powered by either 2 alkaline 1.5V batteries or a 5V Micro USB power adapter.

Cables, power adapters, and other external equipment are not included.

## Connections

There are four connections in the product.

1. External relay outputs and external power supply
2. Not used.
3. External vibrator.
4. Charging connector, Micro-USB.



External output/power supply. Not used External vibrator. Micro-USB charger.

1. Minus
- 2.
3. 8-15V @ 100mA
4. Relay output 4 (NO)
5. Relay output 3 (NO)
6. Relay output 2 (NO)
7. Relay output 1 (NO)
8. Relay common (COM)

## Commissioning

Before usage, 2 AA batteries must be inserted into the battery holder.

Remove the battery cover and insert the batteries as shown.

**Batteries in the set:**



When the product is switched on for the first time, the desired language must be selected.

If the product is switched off, press **Ent** briefly and the product switches on.



Press to move the green arrow to the desired language.

Press **Ent** to select a language.

**NOTE:** After the product has been switched on, when it has been switched off, it may display "Low battery" briefly until it has measured the battery.



## Start screen

Pressing one of the keys turns on the display and shows the time and date.

The user's name, if encoded, is displayed at the top.



### Home screen options:



The name of e.g. a product or a location can be changed to the desired name.



The current time and date are displayed in white text.



The clock must be set.

### Battery status:



Battery OK.



It is recommended that the battery is replaced.



The battery is running low and needs to be replaced.

### Alarm sound:



The alarm sound is switched on.



The alarm sound is switched off.


### Vibrator:



The vibrator is activated.



The vibrator is switched off.

Press  to see the call list.





## Reception of alarms

When the receiver receives an alarm from a coded transmitter, the name, type, and time of the reception and any position are displayed.



Client's name  
Optional text if encoded.  
Type of transmitter.  
Positioning for the alarm.  
Button/function type.  
Relay output 1


The individual texts are only displayed if they are encoded or known.

Press **Ent** to confirm the alarm.

Press **Esc** to reject.

If the alarm call is unanswered, short sounds are heard approximately every five minutes, as notification if not muted.

When a message is received from a transmitter, e.g. with a low battery, the message below is displayed.



Press **Ent** to confirm the alarms.

**NOTE:** "Last call" displays a low battery message for the transmitter where the battery needs to be replaced.



## Emergency call

It is possible to send an emergency call (or alarm) from the recipient to e.g. other users to indicate "I need help" at a client or as an assault alarm.

The emergency call is activated by pressing the KNOP logo.

However, this requires the receiver to be encoded in other receivers.




Emergency call	Emergency call	Emergency call
Transmit Emergency call!	Transmit Emergency call!	Transmit Emergency call!
	OK	Error!
Transmitting	Accepted	Failed



## Batteryalarm

The receiver indicates the battery status on the display and in case of discharged battery with sound indications.

### **Battery status:**

-  Battery OK.
-  Replacing the batteries is recommended.
-  The battery is next used up and needs to be replaced.

If the symbol changes to red, immediately insert the device into the charger.

At the same time, a few short beeps will be heard about every five minutes, indicating low battery, and the screen below will appear:





## Clients and transmitters

Before you can encode a transmitter in the receiver, you must create a client to which the transmitter belongs.

A client can be a person, a room, or something else. For example, the client can be **Sven** (used as an example).

Our examples our client is called **Sven** can have one or more transmitters associated with him.

For example, **Sven** (our client) has a wristband transmitter, a bed alarm, a motion sensor, and other products.

When **Sven** activates an alarm from his products, the display shows that the alarm comes from **Sven** and which product is the sender, so that the cause of the alarm is known from the beginning.

This of course assumes that names etc. are encoded in the desired way.

### Clients created under Settings / Client / New client:



Main menu	Settings	Clients
↑ Hibernate Select Group Vibrator: On Sound: Low → Settings	Ver: 3.8.1 ID: 123456 Clients → Edit group ↓ Set clock	→ Add client <b>Sven</b> All TX's.

### Add under Setting / Client / Client (Sven) / Add TX unit:


Main menu	Settings	Clients	Sven
↑ Hibernate Select Group Vibrator: On Sound: Low → Settings	Ver: 3.8.1 ID: 123456 Clients → Edit group ↓ Set clock	Add client → <b>Sven</b> All TX's.	↑ Add TX unit Select group Client name → Client text Delete



## Main menu

If you press  and  at the same time, either when the screen is off or on the home screen, the main menu opens.

The green arrow on the left shows which line in the menu can be selected.



Press  to move the green arrow either up or to the right in the menu.

Press  to move the green arrow either down or left in the menu.

Press  to leave the menu.

Press  to select the menu item with the green arrow next to it.

**NOTE:** After about 60 seconds of inactivity, the menu closes automatically.


### Latest alarms:




  Scroll between alarms.

**GREEN** Means that the alarm was accepted.  
**RED** Means that the alarm was rejected.  
**YELLOW** Means the alarm went unanswered.

### See more information about the alarms:



 See the client's position and which transmitter the alarm is coming from.


  Scroll through multiple alarms.

When an alarm is received from a transmitter that has, for example, a low battery, the message is displayed.

  to scroll through other alarm calls.



## Delete the call list:



Press **Ent** to confirm.


## Suspend/Hibernate receiver:



**Ent** Confirm or **Esc** to reject.

The receiver is switched on by pressing **Ent**.

## Group selection:



**Ent** Select or deselect a group.


**Esc** Return to menu.

**↓** **↑** Moves arrow to the group you want to change.

**1** Group selected

**1** Group deselected


## Vibrator function:



**Ent** Activate or deactivate the vibrator function.

**Esc** Return to menu.

## Sound setting:



**Ent** Choose between **High**, **Low** or **Off** sound when receiving alarm calls.



## Settings and encoding:

Main menu	Settings	Settings	Settings
↑ Hibernate Select Group Vibrator: Off Sound: High → Settings	Ver: 3.8.1 ID: 123456 → Clients Edit group ↓ Set clock	↑ Set clock Click: Low Product name → Passcode ↓ Outgoing	↑ Passcode Outgoing Positions Language → Turn off

Ver: 3.3.14

The software version of the product.

ID: 58651

Serial number of the product.

### Add new client:

Settings	Clients
Ver: 3.8.1 ID: 123456 → Clients Edit group ↓ Set clock	→ Add client <b>Sven</b> All TX's.

Select **Add client** or select client and press **Ent**.It will ask for the client's name and optional text.  
Remember to add TX devices.

### Change client:

Settings	Clients	Sven
Ver: 3.8.1 ID: 123456 → Clients Edit group ↓ Set clock	Add client → <b>Sven</b> All TX's.	Room 10 → <b>TX901</b> Add TX unit Select group ↓ Client name

Sven
Room 10
→ <b>TX901</b>


Client name.

Optional text, e.g., name of room.


Product already encoded on the client.



## Add transmitters:



Now activate the desired transmitter.




Accept new transmitter.



If the transmitter is already encoded, the menu below appears, allowing you to move a transmitter between clients.


## Add transmitter to groups:



- Select or deselect a group.
- Return to menu.
- Moves arrow to the group you want to change.

- Group selected
- Group deselected

## Change the client's name:



- Choose letter.
- Select letter.
- Delete last letter.

To save the text select symbol and press .

## Add/change client optional text:




- Choose letter.
- Select letter.
- Delete last letter.

To save the text select symbol and press .






## Delete client:




**Ent** Confirm delete.  
**Esc** Return to menu.

## Change encoded transmitter:



Product name of the transmitter.  
 Transmitter ID.  
 The software version of the transmitter.  
 Battery level in the transmitter (if applies).  
 Selected audio signal when alarm is received.


## Change name of the transmitter:



**↓** **↑** Choose letter.  
**Ent** Select letter.  
**Esc** Delete last letter.

To save the text select symbol **◆** and press **Ent**.

## Change alarm signal:



**↓** **↑** Scroll between audio signals.  
**Ent** Select audio signal.

## Change call text:



This is a TX901 example.  
 It is possible to change the shown text for a button etc.  
 In KNOPtool used defined texts can be made.



## Reset call text:

Sven	Edit text	Call texts
↑ Signal 1 Optional text Alarm signal → Call texts ↓ Output	↑ Extern 2 Extern 3 Extern 4 Low battery → Reset texts	Reset texts  Confirm?  <div> <input checked="" type="checkbox"/> <input type="checkbox"/> </div>

Restore call texts to factory setup.

## Relay outputs setting:

Sven	Output	1. Green
↑ Optional text Alarm signal Call texts → Output ↓ Call: Normal	→ 1. Green 2. Red 3. Yellow 4. Blue ↓ KNOP button	No output → Output: 1 Output: 2 Output: 3 Output: 4

Here is a TX901 example.

By default, all alarms and **Low battery** go to output 1.

The different alarm types or **Low battery** from the transmitter can be selected individually.



## Accepted alarm types:

Sven	Sven
↑ Alarm signal	↑ Alarm signal
Call texts	Call texts
Output	Output
→ Call: Normal	→ Call: All
↓ Timer setup	↓ Timer setup

**Ent** Toggle between **Normal** and **All**.

All accepts both types of alarms Normal and Emergency calls.

## Transmitter activation timer:



Sven	Sven
↑ Call texts	→ Timer: Off
Output	Start 1 00:00
Call: All	Stop 1 00:00
→ Timer setup	Start 2 00:00
Delete	Stop 2 00:00

**Ent** Toggle between **On** and **Off**.

**↓** **↑** Change hour/minutes.

**Ent** Change selection.

## Remove unit from client:

Sven	Sven
↑ Call texts	Delete:
Output	TX901
Call: All	Confirm?
Timer setup	 
→ Delete	

**Ent** Delete confirm.

**Esc** Return to menu.



## Group settings:

Settings	Setup group
Ver: 3.8.1	Group: 1
ID: 123456	
Clients	
→ Edit group	↓
↓ Set clock	1 2 3 4 5 6 7 8

- ↓ ↑ Move between groups.
- Ent Select group.

## Group activation timer on/off:

Group: 1	Group: 1
→ Timer: Off	→ Timer: On
Group name	Group name
Timer setup	Timer setup

- Ent Toggle between **On** and **Off**.
- Esc Return to menu.

Don't forget to set the time in "Settings - Set clock".

## Change group name:

Group: 1	Group: 1
Timer: On	Save?
→ Group name	◆ 4 5 6 7 8 9 a b c d
Timer setup	Gruppe: 1
	Enter
	group name

- ↓ ↑ Choose letter.
- Ent Select letter.
- Esc Delete last letter.

To save the text select symbol ◆ and press Ent.

## Change the timer:

Group: 1	Group: 1
Timer: On	Start 07:00
Group name	Stop 15:00
→ Timer setup	

- ↓ ↑ Change hours or minutes.
- Ent Select items.
- Esc Return to menu.

## Set clock:

Settings	Set clock
↑ ID: 123456	2022-10-27
Clients	11:21
Edit group	Automatic
→ Set clock	Year: 2022
↓ Click: Low	

- ↓ ↑ Change item.
- Ent Scroll between **Year**, **Months**, **Time** and **Date**.

**Note:** If the receiver receives a clock update from a KNOP product or USB901, the clock is automatically set accordingly.



## Button click sound:

Settings	Settings
↑ Clients	↑ Clients
Edit group	Edit group
Set clock	Set clock
→ Click: Low	→ Click: High
↓ Product name	↓ Product name

**Ent** Toggle between **Low**, **High** and **Off**.

**Esc** Return to menu

## Change the product (user) name:

Settings	Name
↑ Edit group	Save?
Set clock	◆ 456789 a b c d
Click: High	RX901C
→ Product name	Enter
↓ Passcode	product name

**↓** **↑** Choose letter.

**Ent** Select letter.

**Esc** Delete last letter.

To save the text select symbol ◆ and press **Ent**.

## Change passcode:

Settings	Passcode
↑ Set clock	New code
Click: Low	
Product name	
→ Passcode	Code: 1234
↓ Outgoing	

**↑** **↓** Select digits (0 to 9).

**Esc** Delete last digit.

The code must always be 4 digits.

## Function of the KNOP button:

Settings	Outgoing
↑ Click: Low	→ Call: All
Product name	
Passcode	
→ Outgoing	
↓ Positions	

**Ent** Change function, choose between **All**, **Normal** or **Off**.


## Positions:

Settings	Positions
↑ Product name	Door
Passcode	→ Show position
Outgoing	Add new
→ Positions	Door
↓ Language	
	Door
	Door

Current position.  
Already encoded.




## View position:



Shows the current position of where you are.

## Add position:




**Ent** Save new position

**Esc** Return to menu.

Note that the receiver must be within range of a position field.

## Change position name:




**↓** **↑** Choose letter.

**Ent** Select letter.

**Esc** Delete last letter.

To save the text select symbol ♦ and press **Ent**.

## Remove position:



**Ent** Delete position.

**Esc** Return to menu.

## Languages:




**↓** **↑** choose language.

**Ent** Accept the chosen language.

**Esc** Return to menu.



Turn off the product:



**Ent** Confirm **Turn off**.

**Esc** Return to menu.

Turn the receiver back on, press **Ent**.

**Note:** The clock and date are reset.



## ***Maintenance***

### **Update of RX901C**

It is possible to update the software of the RX901C, as well as to download and transfer settings with the KNOPtool. For further descriptions, please refer to the KNOPtool manual.

## ***Cleaning***

The product can be cleaned with a moist cloth or a wet wipe.

## ***Spare parts and accessories***

Spare parts and accessories can be ordered from Knop Elektronik A/S at [www.knop.dk](http://www.knop.dk).





## Technical data

<b>Frequency HF:</b>	869.2125 MHz
<b>Range:</b>	Up to 1500 m between 2 RX901B <sup>(1)</sup>
<b>Battery:</b>	2 pcs. AA alkaline batteries
<b>Battery time:</b>	Approximately 180 days with 10 alarms per day (use of vibrator can reduce this).
<b>Charge:</b>	None.
<b>Display:</b>	Graphic LCD.
<b>Menu texts:</b>	3 languages: Danish, Norwegian, and Swedish.
<b>Log function:</b>	Saves the last 50 alarm calls.
<b>Memory:</b>	250 transmitters 150 position rod POS901
<b>Surrounding environment:</b>	Indoor use. $\leq 90\%$ non-condensing
<b>Ambient temperature:</b>	0°C to +40°C
<b>Type of cabinet:</b>	Colour RAL9001 ABS
<b>Cabinet dimensions:</b>	W: 65mm, H: 120mm, D: 20mm.
<b>Density:</b>	IP20
<b>Weight:</b>	100g















The right to make changes is reserved.

All rights reserved.

© **KNOP** Elektronik A/S

<sup>(1)</sup> Measured outdoors with an unobstructed view between transmitter and receiver. In buildings, the range is reduced.

## Explanation of symbols and approvals

	<b>This product complies with:</b> Directive 2017/745/EU      MDR Directive 1907/2006/EU      REACH Directive 2011/65/EU      RoHS Directive 2012/19/EU      WEEE ISO 14971:2019      Risk Management for Medical devices EN 301 498-1 V2.2.3      Electro Magnetic Compatibility EN 301 489-3 V3.1.1      Electro Magnetic Compatibility EN 50130-4:2011 + 2014      Immunity alarm systems EN 300 220-1 V3.1.1      Short Range Devices EN 300 220-2 V3.1.1      Short Range Devices EN 300 220-3 V2.1.1      Short Range Devices EN 60601-1-2:2014 + 2015      Electromagnetic compatibility (Medical) EN 62368-1:2020      Electrical safety EN 50581:2012      Hazardous substances
	The product must not be used if the packaging is damaged.
	<b>Medical Device</b> Class 1, rule 1
	<b>Manufacturer</b> KNOP Elektronik A/S, Fabriksvej 20, DK-7600 Struer, Denmark
	Read the manual(s) before installation and commissioning at <a href="http://www.knop.dk">www.knop.dk</a> .
	Interference may occur in an environment with equipment marked with this symbol.
	Must be protected against liquids.
	0 °C to +40 °C, temperature limit for transport/storage and use.
	The product must not be disposed of with normal household waste.
	<b>Single registration number</b> DK-MF-000025631
	<b>Unique Device Identifier</b> 05744002850489
	<b>Product reference</b> RX901C, Receiver with relay outputs
	<b>Serial number</b> Placed on the product
	<b>Ingress Protection Code</b> IP20